CALIBRATION CERTIFICATE

SERIAL NUMBER: DF25A-5863

PAGE 1 **OF** 3

CALIBRATION LABORATORY: Labsphere, Inc., Reflectance Calibration Laboratory

CALIBRATED STANDARD OR SAMPLE

Spectralon Reflectance Target, Sample I. D: SRT-20-050

CALIBRATION REQUESTED

8°/Hemispherical Spectral Reflectance Calibration from 250-2500 nm reported in 50 nm intervals.

APPLICABLE DOCUMENTS

Labsphere #QP-14012-000 General Reflectance Calibration Requirements

#QP-14069-000 Background Correction for RSA-PE-19

#OP-14068-000 Procedure For Calibration of Spectralon Material

#QP-14074-000 Morning Check for RSA-PE-19

TEST AND TEST METHOD

The 8°/Hemispherical Spectral Reflectance was measured for the sample listed above. The calibration was performed with standards traceable to the U.S. National Institute of Standards and Technology (NIST). The reflectance was determined by using a Perkin-Elmer 900 double beam ratio recording spectrophotometer (serial # 70038) equipped with a Labsphere RSA-PELA-1000 integrating sphere accessory (serial # 08220625262). Incident flux at 8° from normal was reflected from the above sample then radiated onto the internal surface of the integrating sphere. The sphere's internal surface radiance was proportional to that of a laboratory working standard used to determine the spectral reflectance.

NATIONAL LABORATORY TRACEABLE STANDARD

SRS-99-020-WS-3 (NIST)

CALIBRATION RESULTS MEETS SPECIFICATIONS

See Table I "8°/Hemispherical Spectral Reflectance Factor" See Graph I "8°/Hemispherical Spectral Reflectance Factor"

LABORATORY ENVIRONMENT

Temperature: $23 \pm 3^{\circ}$ C.

ESTIMATE CALIBRATION UNCERTAINTY

The random uncertainty, expressed by standard deviation, of reflectance measurements performed by Labsphere, Inc. is less than 0.005 for the spectral range 300-2200 nm, and is less than or equal to 0.020 over the spectral range 250-2500 nm.

CALIBRATION INTERVAL

If the product is handled with care, recalibration is recommended after one year of use, at the discretion of the user.

This certificate shall not be reproduced except in full, without the written approval of Labsphere, Inc.

Calibrated by:

Title: Optical Calibration Technician

Date:



P.O. Box 70 No. Sutton, N.H. 03260 Tel 603-927-4266 Fax 603-927-4694

Table I 8°/Hemispherical Spectral Reflectance Factor for SRT-20-050

Wavelength	Reflectance
<u>(nm)</u>	<u>Factor</u>
250	0.210
300	0.201
350	0.199
400	0.201
450	0.204
500	0.209
550	0.214
600	0.219
650	0.224
700	0.228
750	0.234
800	0.239
850	0.245
900	0.255
950	0.259
1000	0.262
1050	0.266
1100	0.271
1150	0.276
1200	0.280
1250	0.285
1300	0.289
1350	0.292


Wavelength	Reflectance
<u>(nm)</u>	<u>Factor</u>
1400	0.298
1450	0.302
1500	0.306
1550	0.310
1600	0.314
1650	0.318
1700	0.322
1750	0.326
1800	0.329
1850	0.335
1900	0.338
1950	0.344
2000	0.346
2050	0.350
2100	0.354
2150	0.359
2200	0.359
2250	0.362
2300	0.363
2350	0.368
2400	0.369
2450	0.373
2500	0.373



P.O. Box 70 No. Sutton, N.H. 03260 Tel 603-927-4266 Fax 603-927-4694



PAGE 3 **OF** 3





P.O. Box 70 No. Sutton, N.H. 03260 Tel 603-927-4266 Fax 603-927-4694

